

SalonSERVE®



Microcurrent Machine

SEQ641



- SELECT - Mode Selection:  
Complete with 2 modes:



Hydration / Interval Mode: In this mode, the electrical currents come in pulses with breaks in between. This helps reduce puffiness by helping the body get rid of excess fluid. It also stimulates the skin cells to produce collagen. So, while it helps with hydration like continuity mode, it's also great for keeping the skin firm and youthful over time.



Collagen / Continuity Mode: This mode delivers a steady stream of electrical currents to the skin. It helps improve blood flow and allows the skin to hold onto moisture better, which can make it look more hydrated. Over time, this mode can also help the skin produce more collagen, which keeps it firm and youthful.



# Introduction to Microcurrent technology

The latest facial toning machines use “microcurrent” technology. Using micro current impulses, they achieve almost instant visible results to tighten, tone, and firm the skin. In fact, the results from microcurrent facial toning can be so dramatic that they are often referred to as “non surgical face lifts”. Originally this type of muscle stimulation was used as a sports physiotherapy treatment but has now been successfully adapted to use its toning and firming effect to recontour and rejuvenate facial muscles.

## How does Microcurrent work?

Microcurrent facial toning machines relay safe, painless impulses to the facial muscles and help overly-tightened muscles to relax and under worked muscles to regain strength. Microcurrent facial toning “re-educates” the facial muscles and helps restore them back to their original tone.

The microcurrent treatment utilises micro current technology which uses gentle electrical stimulation to trigger the body’s natural skin enhancement chemicals at a cellular level. Clinical studies have shown that after twenty days of treatment collagen production can increase by 14%, elastin can increase by 48%, and circulation increases by 38%. Microcurrent treatments can also trigger the body’s production of amino acids and ATP, accelerating cell repair and promotes healthier cell production.

## Client suitability and treatment time

All ages of clients can benefit from this treatment. Younger clients can have microcurrent facials as preventative treatments whilst mature clients will benefit from the rejuvenating effects of the microcurrent.

Microcurrent is suitable for both male and female clients. The treatment will last approximately 1 hour and clients should be prescribed a course of 8 - 12 treatments to start with to achieve the desired results.

Follow up treatments are recommended every 4 - 6 weeks to maintain results.

## Contraindications specific to Microcurrent treatments

- Pacemaker
- Pregnancy
- Epilepsy - with doctors approval

- Diabetes - controlled diabetes is safe to treat
- Contagious skin disease
- Lack of skin sensation or allergies - skin sensitivity test should be undertaken prior to treatment
- Metal pins/plates and/or excessive dental fillings or bridge work
- Recent scar tissue in area to be treated
- Retin A and topical steroids
- Roacutane
- Dermal fillers or collagen treatments. Need to wait at least 2 - 3 weeks after these have been administered.

## The Salonserve Bio Ultimate Facial

The Salonserve Bio Ultimate facial is a non-surgical non-invasive cosmetic technique using micro currents to deliver a low level of electrical impulses into the skin which stimulates chemical reactions within the skin’s cellular level. The Bio Ultimate electrical current replicates the body’s own natural electrical current and provides a number of benefits:

- Improving muscle tone by tightening and firming the sagging muscles
- Refining fine lines
- Defining facial features
- Improving blood and lymph circulation to the face and neck
- Improving the elimination of toxins and waste from the skin
- Stimulating sluggish circulation.
- Improving skin colour
- Increasing cellular regeneration.

## Electrodes used with the Salonserve Bio Ultimate Facial

There are several different types of electrodes that can be used within the microcurrent treatment:

### ■ Facial Bio Rods

Use with disposable paper stem cotton buds. Cut the stem to fit into the end of the rod so just the bud is protruding.

### ■ Body Bio Rods

These are the large rods used for lifting main body muscles.

## The Salonserve Bio Ultimate Treatment Procedure

Work on cleansed and toned skin.

Select the electrodes appropriate for the treatment to be carried out. If using facial rods cut paper stem cotton buds and insert into the ends.

Plug electrode into the socket on the main frame of the machine, as indicated on the faceplate.

Switch the machine on.

Press the On/Off key to activate the programmes. Select the programme according to the desired results and clients skin conditions.

**Body** - body treatments

**Facial** - micro current facial

**Bio** - for clients wanting a more natural lift

Adjust the treatment time by increasing or decreasing the timer keys.

Adjust the current output by increasing or decreasing the Output key.

To begin the treatment start on a low output and work your way up. As the client becomes more accustomed to the treatment the quicker the treatment output can be increased.

Unlike Faradic treatments the micro current shows little or no muscle contraction so ensure that client comfort is maintained throughout the treatment.

Select a suitable conducting gel for the client's skin type and apply an even layer over the face and neck with fingertips or a mask brush. (no need to use gel if using the magic gloves)

■ On second visit increase output

■ On third visit increase output again

■ On fourth and subsequent visits increase to maximum output

Always check the client is comfortable with the increased output levels.

There may be evidence of erythema but this is perfectly normal as it is a sign of increased circulation.

Work slowly and rhythmically covering all the areas thoroughly.

On completion of the treatment remove the electrodes and switch the machine to off.


Remove any remaining conductive gel using damp cotton wool or warm water and sponges.

To complete the treatment tone and moisturise the skin.

# Notes

A series of horizontal dotted lines for taking notes.

1



Start with rods in the centre of the muscles and slide them along the muscle in opposite directions. When placing the rods press firmly then glide away from each other. Move across the neck avoiding the windpipe.

2




Work either side of the windpipe sliding the rods out towards the Sternocleidomastoid muscles.

3




Zig zag rods upwards towards the jaw line avoiding the windpipe.

4



Work along the jaw line gliding the rods towards the ear.

5



Zig zag rods along the jaw line towards the ears.

6




Start with the one rod under the eye and the other by the mouth. Press the rods to the skin and then glide towards each other pressing the muscles between the rods. Hold the muscle then release.

7



Zig zag rods across the cheeks towards the ears.

8




As on the neck use the rods to stretch the corrugator muscle. Place the rods firmly in the centre of the muscle and then glide away from each other.

9



The same application is carried out across the frontalis muscle. Start in the centre of the forehead and work outwards.

10




Zig zag rods across the forehead.

11



Place the rods at the side of the eye and glide away from each other horizontally and vertically.


12




Zig zag rods around the orbicularis oculi.

Movements should be slow and rhythmical with a firm pressure. Always check pressure with client.


# Sequence for Microcurrent application on the breast




Apply the lifting gel over the chest that joins the top of the breast (pectoralis muscles) across to the underarm (axilla). Apply the gel below the breast where the tissues join the top of the abdominal wall (serratus anterior muscle).




Work from the top of the breast above the nipple to just below the clavical. Use slow rhythmical vertical lines from axilla to breast bone but do not work on the breast bone. Start with Body Rods together and smooth away from each other.



Work on the top of the breast using horizontal movements as above from breast bone to axilla.




Work just below the breast area. Use slow, rhythmical vertical and horizontal movements from breast bone to axilla. On larger breasted clients you may need to ask that the breast is lifted for this application.




As No. 4 above.




# Sequence for Microcurrent application on the buttocks




Apply the lifting gel over the buttocks (gluteal muscles) and from the middle of the back of thighs (hamstring muscles) to under gluteal ridge.




Work from middle of buttock out towards lower back and gluteal ridge. With the body rods use firm vertical stretching movements.



Using circular movements with the rods work from middle of buttocks out towards outer buttock and gluteal fold.




Work just below the gluteal ridge on the upper back of thighs (hamstrings). Using body rods start from mid thigh and gluteal ridge and bring rods together with firm vertical movements.




Starting with rods under the gluteal ridge at the top of upper thigh firmly work from centre to outer and inner thigh.


# Sequence for Microcurrent application on the abdomen




Apply the lifting gel over the abdomen (rectus abdominus and transversus muscles) and from the centre of the abdomen to under the rib cage and from centre of abdomen to waist.




Work from middle of abdomen out towards pubic area and under rib cage. With the body rods use firm vertical stretching movements.



Starting with the rods together at the centre of the abdomen use horizontal stretching movements out towards the waist and side of abdominal area.




Using circular movements with the rods work from middle of abdomen out towards the waist and side of abdominal area.



Starting with rods at the centre of the abdomen and pubic area bring them together lifting and pressing lower rectus abdominal muscles together.


# Sequence for Microcurrent application on the thighs

1




Apply the lifting gel all over the thighs from the top of the knee (patella) to the groin.

2



Work from middle of thighs out towards groin and above patella. With the body rods use firm vertical stretching movements.

3



Using circular movements with the rods work from middle of thighs out towards the groin and above patella.

## Technical Specifications

Power supply: 220 volt  
Hertz: 50Hz  
Watts: 15W

## Machine Care

- Do not immerse any parts of the machine or handset in water.
- Clean the electrodes with surgical spirit on cotton wool between each client.
- Clean the unit with a damp soft cloth.

SalonServe Hair & Beauty Supplies Ltd, Lexicon  
House, Station Road, Irthlingborough, NN9 5QF

Telephone: 01933 410022

[www.salonservice.com](http://www.salonservice.com)